# ARTIA™ Reconstructive Tissue Matrix

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## DEVICE DESCRIPTION

ARTIA™ Reconstructive Tissue Matrix (ARTIA™ Tissue Matrix) is a surgical mesh that is derived from porcine skin and is processed and preserved in a patented aqueous phosphate buffered solution containing matrix stabilizers. ARTIA™ Tissue Matrix is intended to perform as a surgical mesh for soft tissue repair while presenting a scaffold to the patient's tissue. ARTIA™ Tissue Matrix consists of a sterilized sheet of processed porcine dermis provided in various geometric configurations and packaged in a plastic holder enclosed within a pouch.

Use of the surgical mesh provides for a strong and biocompatible implant and will incorporate into the patient's tissue with associated cell and microvascular ingrowth.

## INDICATIONS

ARTIA™ Tissue Matrix is intended for use as a soft tissue patch to reinforce soft tissue where weakness exists and for the surgical repair of damaged or ruptured soft tissue membranes. Indications for use include repairs in breast plastic and reconstructive applications which require the use of reinforcing material to obtain the desired surgical outcome.

It is indicated to be used to reconstruct, to recontour and to reform the host's human soft connective tissue particularly where loss of tissue has occurred and as a supporting tissue in breast plastic and reconstructive surgical applications.

ARTIA™ Tissue Matrix is supplied sterile and is intended for single patient one-time use.

## CONTRAINDICATIONS

- This surgical mesh is derived from a porcine source and should not be used in patients with known sensitivity to porcine material.
- Polysorbate 20 is a component of the aqueous phosphate buffered solution and therefore ARTIA™ Tissue Matrix should not be used in patients with a known sensitivity to this material.

## WARNINGS

- **Do not resterilize.**
- Do not use if the package is opened or damaged. Do not use if seal is broken or compromised. Do not use if the temperature monitoring device does not display "OK".
- After use, handle and dispose of all unused product and packaging in accordance with accepted medical practices and applicable national and regional environmental laws on disposal of packaging and biological materials.
- ARTIA™ Tissue Matrix cannot be reused once it has been removed from the packaging and/or is in contact with a patient without increased risk of patient-to-patient contamination and subsequent infection.

## PRECAUTIONS

- Discard surgical mesh if handling has caused possible damage or contamination.
- Discard surgical mesh if it is past the use-by-date of the product (indicated as 4 digit year, 2 digit month, and 2 digit day [YYYY-MM-DD]).

## ARTIA™ Reconstructive Tissue Matrix

### Instructions For Use | English

- Ensure that the surgical mesh is put into a sterile basin and covered with room temperature sterile saline or room temperature sterile lactated Ringer's solution for a minimum of 2 minutes prior to implantation.
- Place surgical mesh in maximum possible contact with healthy, well-vascularized tissue to promote cell ingrowth and tissue remodeling.
- ARTIA™ Tissue Matrix should be hydrated and moist when the package is opened. If ARTIA™ Tissue Matrix is dry, do not use.
- If a tissue punch-out piece is visible, remove using aseptic technique before implantation.
- Presence of a significant microbial load may affect overall performance of surgical mesh. Utilize bioburden-reducing techniques to minimize contamination levels at the surgical site, including, but not limited to, appropriate drainage, debridement, negative pressure therapy, and/or antimicrobial therapy.
- Certain considerations should be made in order to reduce the risk of adverse events when performing breast reconstruction using a surgical mesh such as ARTIA™ Tissue Matrix. Please see the following sections for more information: Product/Patient Selection, Technique Guidance, and Post-Operative Care.

## ADVERSE EVENTS

Potential adverse events are those typically associated with the surgical mesh materials and/or their implantation procedures including, but not limited to, infection, foreign body response, failure to integrate, hematoma, seroma formation, inadequate tissue perfusion (including tissue ischemia and incision or tissue necrosis), wound dehiscence, lack of support to host tissue, and inflammation. If an unanticipated event occurs, alteration of surgical plan may be necessary at the surgeon's discretion.

## STORAGE

- ARTIA™ Tissue Matrix is a sterile medical device that should be stored in a clean, dry location at -8°C to 30°C.
- Refer to the temperature monitor located on the product carton to ensure that the product has been stored within its temperature limits. Only use the product if the included temperature monitor displays "OK" on the screen. If the screen displays anything other than "OK", do not use the product.
- It is to be stored in its original packaging.
- The use-by-date of the product is indicated as 4 digit year, 2 digit month, and 2 digit day (YYYY-MM-DD).

## STERILIZATION

This product is sterilized by E-Beam.

## PATIENT/PRODUCT SELECTION

The following considerations are intended to serve only as general guidelines. They are not intended to supersede institutional protocols or professional clinical judgment concerning patient care. Physicians should be familiar with surgical procedures and techniques involving biological surgical mesh before using ARTIA™ Tissue Matrix.

- Carefully consider the risk/benefit balance of using a surgical mesh such as ARTIA™ Tissue Matrix in patients with significant comorbidities; including but not limited to, obesity, smoking, diabetes, immunosuppression, malnourishment, poor tissue oxygenation (such as COPD), and pre- or post-operative radiation.
- Carefully consider patient factors (e.g. skin flap size and viability) and product sizing when selecting a surgical mesh such as ARTIA™ Tissue Matrix. It is important to create space for the tissue expander/breast implant that precisely fits the tissue expander/breast implant. This will help to avoid any potential space for fluid accumulation within the implant pocket. With proper product sizing the skin flaps should precisely fit over the muscle/ARTIA™ Tissue Matrix layer for sub-muscular coverage or ARTIA™ Tissue Matrix layer for pre-pectoral coverage without redundancy and without excessive tension.

## INSTRUCTIONS FOR PREPARING ARTIA™ TISSUE MATRIX FOR SURGICAL USE

These instructions are designed to serve only as a general guideline. They are not intended to supersede institutional protocols or professional clinical judgment concerning patient care. Physicians should be familiar with surgical procedures and techniques involving biological surgical mesh before using ARTIA™ Tissue Matrix.

## REQUIRED MATERIALS

- Sterile forceps.
- Soaking fluid: room temperature sterile saline or room temperature sterile lactated Ringer's solution.
- One sterile basin per piece of ARTIA™ Tissue Matrix.

## PREPARATION

1. Open the carton and remove the pouch.
2. Peel open the pouch and remove the plastic holder using aseptic technique. The plastic holder is sterile and may be placed directly into the sterile field.
3. Open the plastic holder carefully and aseptically remove the surgical mesh. Always use sterile gloves or forceps when handling ARTIA™ Tissue Matrix.
4. Place the surgical mesh in a sterile basin and submerge in room temperature sterile saline or room temperature sterile lactated Ringer's solution for a minimum of 2 minutes.
5. Keep the surgical mesh in the room temperature sterile solution until ready for implantation. Surgical mesh can be retained in sterile solution for a maximum of 4 hours.

## IMPLANTATION

1. Prepare the surgical site using standard techniques.
2. ARTIA™ Tissue Matrix may be folded, trimmed or cut as required to fit the surgical site using aseptic technique, ensuring allowance for overlap.
3. Transfer ARTIA™ Tissue Matrix to the surgical site using sterile gloves or forceps.
4. Suture ARTIA™ Tissue Matrix into place.
5. Complete the standard surgical procedure.
6. Discard any unused portions of ARTIA™ Tissue Matrix as per institutional procedures.

The following suggestions are not intended to supersede institutional protocols or professional clinical judgment concerning patient care. Physicians should be familiar with surgical procedures and techniques involving biological surgical mesh before using ARTIA™ Tissue Matrix.

## TECHNIQUE GUIDANCE


- As with any surgical implant, careful aseptic technique should be practiced and contact of the mesh with the patient's skin should be minimized.
- The quality and adequate thickness of the skin flap should be carefully assessed to determine the best method of breast reconstruction for the patient. Carefully assess the skin flaps to ensure that they are well perfused prior to proceeding with breast reconstruction using ARTIA™ Tissue Matrix.
- ARTIA™ Tissue Matrix should be covered with healthy, viable tissue; portions of the skin flap that appear ischemic or necrotic should be addressed using clinical judgment.
- Care should be taken to avoid excessive tension on the skin flaps at the time of closure.
- The subcutaneous space created for the tissue expander/breast implant should precisely fit the tissue expander/breast implant so as to avoid any potential space for fluid accumulation within the implant pocket. For sub-muscular coverage, with proper product sizing the skin flaps should precisely fit over the muscle/ ARTIA™ Tissue Matrix layer without redundancy and without excessive tension. For pre-pectoral coverage, with proper product sizing the skin flaps should precisely fit over the ARTIA™ Tissue Matrix layer without redundancy and without excessive tension. Additional caution is recommended for larger implants (>500ml) as the weight and size of the implant may contribute to the tension on the skin closure.
- It is recommended that drains, if used, should be placed in the tissue expander/breast implant space as well as in the space between the skin flap and the ARTIA™ Tissue Matrix so as to reduce fluid accumulation.

## POST-OPERATIVE CARE

- Proper post-operative care may include the use of appropriate drainage and negative pressure therapy.

- Drains should be left in place until the amount of drainage reaches a level of less than 30ml per 24 hour period for 3 consecutive days. Typical timeline for drain removal is generally 1 to 2 weeks, but should be customized to the clinical circumstance.
- As with any post-operative care, aseptic technique should be practiced when required to minimize contamination to the surgical wound.

## DEFINITIONS

	CONTENTS		OPEN HERE		DO NOT USE IF PACKAGE IS DAMAGED
	CORRECT SETTING		INCORRECT SETTING		MANUFACTURED BY
	ATTENTION, SEE INSTRUCTIONS FOR USE		STERILIZED USING IRRADIATION		DO NOT RE-STERILIZE
	FOR SINGLE USE ONLY		TEMPERATURE LIMITATIONS		USE BY DATE
	CATALOG NUMBER		LOT NUMBER		AUTHORIZED REPRESENTATIVE IN THE EUROPEAN COMMUNITY
	OVOIDAL		CONTOUR		RECTANGLE
	PERFORATED	<b>RxONLY</b>		CAUTION: FEDERAL (USA) LAW RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN	

## CUSTOMER SUPPORT

For product complaints and potential adverse events, please contact your local Sales Representative, or +44-1628-497-456.

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Allergan® and its design are trademarks of Allergan, Inc.
ARTIA™ is a trademark of LifeCell Corporation, an Allergan affiliate.
© 2021 Allergan. All rights reserved.
LifeCell Corporation
Allergan Pharmaceuticals International Limited